

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicaids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company



EP 90

Electric cum Manual Operated Suction Unit
A must for places where electric supply is uncertain.

CE
0434

Housing: MS Powder Coated with top and panel of rust proof and corrosion resistant moulded ABS.

Capacity: Electric -700 ± 10mm Hg, manual -600 ± 10mm Hg with a per stroke average volume of 200ml.

Pump Type: Rotary vane pump (Elec.) and Piston (manual)

Jars: Autoclavable PC jars (2x2 Ltr.), 250 x 55mm tray and a catheter holder with two size hooks to hold patient tube & catheters.

Vacuum Gauge: 6.25 cms. 0-760mm Hg calibration.

Power: 230V, 50 Hz, 2 ± 0.5 Amps, 200 watts. (110V on request)

Noise Level: 50 dB A ± 3 at 1.2 mtr.

Dimension & Weight: 40 x 40 x 90 cms: 30 Kg.

Additional safety jar, fuse & 2 volt foot switch.

Anand Medicaids Pvt. Ltd.

16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com