

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicaids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company

EURO-DEL

Ideal Delivery Aid for Obstetricians : Safe for Mother/Baby

CE
0434

Housing: Engineering moulded ABS plastic cabinet with built-in handle and jar hangers

Pump Type: Oil-free Piston Pump

Jar: 1 ltr PC

Capacity: 0-575 mm Hg \pm 10. Fine regulation valve for necessary adjustment with break mechanism under thumb in ventouse cup. Hi/Lo switch provides full or half vacuum with one flick

Vacuum Gauge: Bourden type 5 cms. 0-760mm Hg calibration.

Filter: Bacterial filter

Power: 220 V, 50/60 Hz, Single phase. (110V on request)

Ventouse cup: Silicone - 60 mm and 70 mm : minimizes risk of trauma and cephalohematome of fetal cranium

Dimension & Weight: 33 x 14 x 29 cms. 6 kg



Anand Medicaids Pvt. Ltd.

16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com