

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicoids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company

EUROVAC PLUS

The Portable Powerhouse Pump

CE
0434

Housing: Robust engineering plastic moulded, shock-resistant ABS with 2 size hooks to hold tube.

Capacity: -700mm Hg \pm 10 Vacuum, 35 ltrs/min.

Pump Type: Double diaphragm Pump

Jars: Autoclavable 2 ltr. collection jar (optional 2nd collection jar on the other side), safety jar with mechanical & electronic overflow monitoring system and bacterial filter. Silicone tube 7 (x11) mm x 2 mtr.

Vacuum Gauge: 6.25 cms. 0-760mm Hg (also kPa) calibration

Power: 220 V AC, 50/60 Hz, 125 watts. (110V on request)

Dimension & Weight: 20 x 42 x 36 cms. 8.5 kg.



Anand Medicoids Pvt. Ltd.

16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com