

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicaids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company



Eurovac-Elite

Anand's Oil Free Diaphragm Pump on Trolley

CE
0434

Housing: Engineering moulded ABS plastic box

Capacity: -700mm Hg \pm 10 Vacuum, 35 ltrs/min.

Pump Type: Double diaphragm pump

Jars: 2 x 2.5 Ltr jars with two-way changeover lever and safety jar with mechanical and electronic monitoring system and bacterial filter

Vacuum Gauge: 6.25 cms (2.5 inch, square dial). 0-760mm Hg (also kPa) calibration.

Power: 220 V AC, 50/60 Hz, 125 watts. (110V on request)

Dimension & Weight: 39 x 55 x 119 cms. 18 kg (with trolley)

Adaptable to hold Swiss Medela Autoclavable safety jar with easy filter replacement.

Controls on working height with extra pedal-operatable switch.

Trolley with non-marking castors having brake facility. Choice of 1 Ltr. jars and of using them singly or in tandem.

Anand Medicaids Pvt. Ltd.

16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com