

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicoids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company



HI-VAC

**Anand's High Vacuum Suction Unit - India's
First ISI Marked Suction Apparatus.**

IS:4533



CM/L-9180073

CE
0434

Housing: Spot Welded, Oven baked cabinet

Capacity: -700 mm Hg \pm 10 regulable, flutter free vacuum control knob, 45 Ltrs. / min.

Pump Type: Double rotary vane pump.

Jars: Wide mouthed 2 x 2 Ltrs. Self sealing bungs with overflow safety device.

Vacuum Gauge: Bourden type 10cm. Dia. 0-760mm Hg (also kg/cm²) calibration.

Power: 230V, 50 Hz, 4 \pm 0.5 Amps, 430 watts. (110V on request)

Noise Level: 50 dB A \pm 3 Almost whispers.

Dimension & Weight: 50 x 38 x 70cms, 40Kg.

#16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com