

WHEN **BETTER UNITS** ARE MADE,

**Anand**<sup>®</sup>

WILL MAKE THEM.

**Anand Medicoids Pvt. Ltd.**

An ISO 13485:2003 & 9001:2000 Certified Company



## HI-VAC JR.

**Anand's High Vacuum Suction Unit - India's  
First ISI Marked Suction Apparatus.**

IS:4533



CM/L-9180073



0434

**Housing:** Spot Welded, Oven baked cabinet

**Capacity:** -700mm Hg  $\pm$ 10 regulable, flutter free vacuum control knob, 28 Ltrs./min.

**Pump Type:** Rotary vane pump.

**Jars:** Wide mouthed 2 x 2 Ltrs. self sealing bungs with overflow safety device.

**Vacuum Gauge:** Bourden type 6.25 cm. Dia. 0-760mm Hg (also inches) calibration.

**Power:** 230 V, 50 Hz, 2  $\pm$  0.5 Amps, 200 watts. (110V on request)

**Noise Level:** 50 dB A  $\pm$  3 Almost whispers.

**Dimension & Weight:** 39 x 33 x 65 cms, 25Kg.

**Anand Medicoids Pvt. Ltd.**

#16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com