

WHEN **BETTER UNITS** ARE MADE,

**Anand**<sup>®</sup> WILL MAKE THEM.

**Anand Medicoids Pvt. Ltd.**

An ISO 13485:2003 & 9001:2000 Certified Company

## INSTAD

Instant Suction Aid - AC/DC - Whenever, Wherever

CE  
0434



**Housing:** Engineering moulded ABS plastic cabinet with built-in handle and jar hangers

**Pump Type:** Oil-free piston pump

**Jar:** 1 ltr PC collection jar with overflow cut-off valve ( 2nd jar optional) Also motor - shut sensor-fitted mini jar

**Capacity:** 0-575 mm Hg  $\pm$ 10 regulable at 12 ltrs/min.

**Vacuum Gauge:** Bourden type 5 cms. 0-760mm Hg calibration.

**Filter:** Bacterial filter

**Power:** 220 V, AC and inbuilt maintenance-free 12 VDC, 2.2 Ah battery, or ambulance /car socket (cable provided). (110V on request)

**Dimension & Weight:** 33 x 14 x 29 cms. 7 kg

**Battery Back up:** 120 minutes on full charge

Supplied in an attractive carry bag

**Anand Medicoids Pvt. Ltd.**

# 16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com