

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicoids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company

OP 108

Anand's Oil - Free Diaphragm Pump

CE
0434



Housing: Engineering moulded ABS plastic cabinet with cavity for catheter holder.

Capacity: -700mm Hg \pm 10 regulable, flutter free vacuum control knob, 32 Ltrs./min.

Pump Type: Double diaphragm pump.

Jars: Autoclavable wide mouthed 2 x 2 Ltrs, self sealing bungs with overflow safety device. Also motor-shut sensor-fitted mini jar and bacterial filter.

Vacuum Gauge: Bourden type 6.25 cm square dial, 0-760mm Hg (also kPa) calibration.

Power: 220 V, 50/60 Hz, 125 watts. (110V on request)

Noise Level: 50 dB

Dimension & Weight: 46 x 34 x 61 cms. 13 kg.

Anand Medicoids Pvt. Ltd.

#16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com