

WHEN **BETTER UNITS** ARE MADE,

Anand[®]

WILL MAKE THEM.

Anand Medicoids Pvt. Ltd.

An ISO 13485:2003 & 9001:2000 Certified Company

P-9

Anand's Non-Tiring Manual-Powered Suction Unit.

CE
0434

Housing: Engineering plastic moulded parts - non corrosive, easy to dismantle, clean, disinfect and reassemble.

Capacity: Outright vacuum -600 mm Hg \pm 10 with average per stroke volume of 200 ml.

Pump Type: Efficient piston pump for creating vacuum instantly.

Jar: Autoclavable 1 ltr. PC jar, silicone patient tube (7 mm ID X 2 mtr).

Vacuum Gauge: Bourden type, 5 cms 0 - 760 mm Hg & 0-100kPa

Dimension & weight: 32 x 17 x 30cms. 2.5 kg.

Portable pedal unit with top view manometer visible to the operator.

Bacterial filter (0.3 micron) with reusable / autoclavable housing between jar and pump.

Ease of operation - one person can operate unaided.

Also useful for field hospitals, midwives, medical centres, mines, aircraft, trains, mountaineers, on-foot paramedics etc.



Anand Medicoids Pvt. Ltd.

16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA

Ph: 91-11-25225225, 25229206, 52464364, Fax: 91-11-25225062

Email: sales@anandind.com, Web: www.anandind.com