



**ANAND MEDICAIDS PRIVATE LIMITED**  
An ISO 13485:2003 & 9001:2000 Certified Company

CE  
0434

## ELECTRIC CUM MANUAL OPERATED SUCTION UNIT ◀



EP 90

**A must for places where electric supply is uncertain.**

**Housing:** MS Powder Coated with top and panel of rust proof and corrosion resistant moulded ABS.

**Capacity:** Electric -700 ± 10 mm Hg, at 25 Ltrs / min., manual - 600 ± 10 mm Hg with a per stroke average volume of 200 ml.

**Pump Type:** Single Rotary vane pump (Elec.) and Piston (Manual)

**Jars:** Wide mouthed autoclavable 2 x 2 Ltr. Jars (Polycarbonate), 250 x 55 mm tray and a catheter holder with two size hooks to hold patient tube and catheters.

**Tubing:** 10 mm ID x 2 mtr (PVC)

**Vacuum Gauge:** 6.25 cms. Dia. 0-760 mm Hg calibration.

**Power:** 230 V, 50 Hz, 2 ± 0.5 Amps, 200 watts. (110 V on request)

**Noise Level:** 50 dB A ± 3 at 1.2 mtr.

**Dimension & Weight:** 40 x 40 x 90 cms. 30 Kg.

Additional safety jar, fuse & 2 volt foot switch.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)  
HI-VAC PLUS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTEUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com