



# ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE  
0434

## AC / DC SUCTION UNITS ◀

### EUROVAC-A



**Anand's Portable AC / DC Suction Unit for Hospitals, Ambulance, Camps, Travel, Homecare etc.**

**Housing:** MS powder coated cabinet

**Capacity:** -575 mm Hg  $\pm$ 10 Vacuum 14 ltrs / min.

**Pump Type:** Oil-free piston pump

**Jar:** Autoclavable 0.5 ltr. collection jar (Polycarbonate) with mechanical / electronic overflow motor shut-off system with autoclavable bacterial filter.

**Tubing:** 8 mm ID X 2 mtr (Silicon)

**Vacuum Gauge:** Bourden type 5 cms. 0-760 mm Hg calibration.

**Power:** Operates on 220 V AC and in-built maintenance-free 12 VDC 4.5 Ah (2 x 6 V, 4.5 Ah) battery, 55 watts. (110V on request)

**Dimension & Weight:** 45 x 15 x 25 cms. 6 Kg.

**Battery Backup:** 120 minutes on full charge

Provided with AC adapter and cigarette lighter socket cable for use in ambulance / car. Can be hung in ambulance or on wall or carried on shoulder in attractive shoulder bag.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)  
HI-VAC PLUS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTEUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407