



ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE
0434

PORTABLE SUCTION UNITS ◀

EUROVAC



Anand's Real Portable Unit with Shock Resistant ABS Case.

Compact, Lightweight, Powerful and Non-Corrosive

Housing: Engineering plastic moulded, shock-resistant ABS with 2 size hooks to hold tube.

Capacity: -575 mm Hg \pm 10 at 25 ltrs / min.

Pump Type: Oil-free diaphragm pump

Jar: Wide mouthed autoclavable 2 ltr. collection jar (Polycarbonate), safety jar with mechanical / electronic overflow motor shut-off system with reusable / autoclavable bacterial filter.

Tubing: 8 mm ID X 2 mtr. (Silicon)

Vacuum Gauge: 6.25 cms (2.5 inch). 0-760 mm Hg calibration.

Power: 220 V AC, 50 / 60 Hz, 125 watts. (110 V on request)

Dimension & Weight: 23 x 42 x 36 cms. 7.5 kg.

Good for different clinical / medical / surgical procedures.

Necessary tools for maintenance / replacement and spares provided.

Product Range

HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)
HI-VAC PLUS (30/65 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC (JR.)
EUROVAC-ELITE
EUROVAC PLUS

MOBILE SUCTION UNIT

MB 36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

ACCESSORIES

POLYCARBONATE JARS
SILICON VENTEUSE CUP
YANKEUR SUCTION TUBE
KARMAN TYPE CANNULAE
BACTERIAL FILTER

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.
New Delhi-110026, INDIA
Ph. +91-11-25225225, 25229206, 42464264
Fax. 91-11-25225062

Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com