



# ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

IS:4533  
CM/L-9180073  
CE  
0434

## HIGH VACUUM SUCTION UNITS ◀

### HI-VAC (JR.)



#### Anand's High Vacuum Suction Unit - India's First ISI Marked Suction Apparatus.

**Housing:** Spot Welded, Oven baked cabinet (Also available with Stainless Steel Body)

**Capacity:** -700 mm Hg  $\pm$ 10 regulable, flutter free vacuum control knob, 28 Ltrs / min.

**Pump Type:** Single rotary vane pump

**Jars:** Wide mouthed 2 x 2 Ltrs. (Polycarbonate) with self sealing bungs and mechanical over flow safety device.

**Tubing:** 10 mm ID x 2 mtr (PVC)

**Vacuum Gauge:** Bourden type 6.25 cm. Dia. 0-760 mm Hg calibration.

**Power:** 230 V, 50 Hz, 4  $\pm$  0.5 Amps, 200 watts. (110 V on request)

**Noise Level:** 50 dB A  $\pm$  3 Almost whispers.

**Dimension & Weight:** 39 x 33 x 65 cms, 25 Kg.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)  
HI-VAC PLUS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTOUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com