



ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE
0434

HIGH VACUUM SUCTION UNITS ◀



HI-VAC PLUSS (30 / 60 LTR / MIN)

Anand's High Vacuum - High Flow Suction Unit

Housing: Plastic moulded cabinet

Capacity: -600 mm Hg \pm 10 mm regulable, flutter free vacuum control knob, 30 / 60 ltrs / min

Pump Type: Double diaphragm pump

Jars: 2 x 2.5 ltr. wide mouthed (Polycarbonate) with changeover lever. Provided with handle for hygiene while emptying / cleaning. Safety jar with mechanical / electronic overflow motor shut off system and autoclavable bacterial filter.

Tubing: 10 mm ID x 2 mtr. (PVC)

Vacuum Gauge: Bourden type 6.25 cm dia 0-760 mm Hg calibration.

Power: 220 V, 50 Hz, Single Phase

Noise Level: 60 dB A \pm 3

Dimension & Weight: 46 x 32 x 86 cms, 18 Kg.

Product Range

HIGH VACUUM SUCTION UNITS

HI-VAC PLUSS (90 LTR/MIN)
HI-VAC PLUSS (30/65 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC (JR.)
EUROVAC-ELITE
EUROVAC PLUS

MOBILE SUCTION UNIT

MB 36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

ACCESSORIES

POLYCARBONATE JARS
SILICON VENTEUSE CUP
YANKEUR SUCTION TUBE
KARMAN TYPE CANNULAE
BACTERIAL FILTER

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.
New Delhi-110026, INDIA
Ph. +91-11-25225225, 25229206, 42464264
Fax. 91-11-25225062

Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com