



# ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE  
0434

## HIGH VACUUM SUCTION UNITS ◀

### HI-VAC PLUSS (90 LTR / MIN)



#### Anand's High Vacuum - High Flow Suction Unit

**Housing:** Plastic moulded cabinet

**Capacity:** -700 mm Hg  $\pm$  10 mm regulable, flutter free vacuum control knob, 90 ltrs / min

**Pump Type:** Double piston pump

**Jars:** 2 x 2.5 ltr. wide mouthed (Polycarbonate) with changeover lever. Provided with handle for hygiene while emptying / cleaning, safety jar with mechanical / electronic overflow motor shut off system and autoclavable bacterial filter.

**Tubing:** 10 mm ID x 2 mtr. (PVC)

**Vacuum Gauge:** Bourden type 6.25 cm dia 0-760 mm Hg calibration.

**Power:** 220 V, 50 Hz, Single Phase

**Noise Level:** 60 dB A $\pm$ 3

**Dimension & Weight:** 46 x 32 x 86 cms, 22 Kg.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUSS (90 LTR/MIN)  
HI-VAC PLUSS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTEUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com