

WHEN **BETTER UNITS** ARE MADE,

**Anand**<sup>®</sup> WILL MAKE THEM.

**Anand Medicaids Pvt. Ltd.**

An ISO 13485:2003 & 9001:2000 Certified Company



## **MB 21**

### **Electric Operated Suction Unit**

CE

**Housing:** Plastic paneled

**Capacity:** -675 Hg (0.09 Mpa) – 25 ltrs/min.

**Pump Type:** Oil-free piston pump.

**Jars:** Wide Mouthed 2x2 ltr autoclavable jars (Polycarbonate) with self sealing bungs and mechanical over flow safety device

**Tubing:** 8 mm ID x 2 mtr (Silicon)

**Power:** 220 V, 50 Hz, 120 watts

**Noise Level:** 50 dB A

**Dimension & Weight:** 38x30x47 cms. 12 Kg.

# 16, Road #33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
Ph: 91-11-25225225, 25229206, 42464264, Fax: 91-11-25225062  
Email: sales@anandind.com, Web: www.anandind.com