



ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE
0434

MOBILE SUCTION UNITS ◀

MB 36



Electric Operated Suction Unit

Housing: Spot Welded, Oven baked cabinet

Capacity: -700 mm Hg \pm 10 regulable, flutter free vacuum control knob, 25 Ltrs / min.

Pump Type: Single rotary vane pump

Jars: Wide mouthed 2 x 2 Ltrs. (Polycarbonate) with self sealing bungs and mechanical over flow safety device.

Tubing: 8 mm ID x 2 mtr (PVC)

Vacuum Gauge: Bourden type 6.5cm. Dia. 0-760 mm Hg calibration.

Power: 230 V, 50 Hz, 2 \pm 0.5 Amps, 200 watts. (110 V on request)

Noise Level: 50 dB A \pm 3 Almost whispers.

Dimension & Weight: 43 x 30 x 68 cms. 27Kg.

Product Range

HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)
HI-VAC PLUS (30/65 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC (JR.)
EUROVAC-ELITE
EUROVAC PLUS

MOBILE SUCTION UNIT

MB 36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

ACCESSORIES

POLYCARBONATE JARS
SILICON VENTEUSE CUP
YANKEUR SUCTION TUBE
KARMAN TYPE CANNULAE
BACTERIAL FILTER

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.
New Delhi-110026, INDIA
Ph. +91-11-25225225, 25229206, 42464264
Fax. 91-11-25225062

Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Ph. 91-1276-267061, 268407

sales@anandind.com, www.anandind.com