



# ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE  
0434

## AC / DC / PEDAL SUCTION UNIT ◀

### MULTIVAC



Anand's Portable AC / DC / Pedal Suction Unit for Ambulances, Hospitals, Camps, Trains, Home care etc.

**Housing:** MS Powder coated cabinet

**Capacity: Electrical:** 0-575 mm Hg  $\pm$  10, Vacuum 14 ltr / min.

**Pedal:** 150 ml per stroke volume

**Pump Type:** Oil free Piston Pump (Electrical) & Piston pump (Pedal)

**Jar:** Wide mouthed autoclavable 2 ltr collection jar (Polycarbonate), safety jar with mechanical / electronic overflow motor shut off system and bacterial filter.

**Tubing:** 8 mm ID x 2 mtr. (Silicon)

**Vacuum Gauge:** Bourden type 5 cm Dia

**Power: Electrical:** Operates on 220 V AC and Pedal (manual). In-built maintenance free 12 V DC battery. 55 Watts.

**Dimension & Weight:** 44 x 28 x 27 cms, 9.5 Kg.

**Battery Backup:** 120 minutes on full charge.

Provided with AC adaptor and cigarette lighter socket cable for use in ambulance / car. Can be hung in ambulance or on wall.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)  
HI-VAC PLUS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTEUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407