



AC / DC / PEDAL SUCTION UNIT



MULTIVAC

- Housing** : MS Powder coated cabinet with electronic overflow protection safety jar.
- Capacity** : Electrical : 0-575 mm Hg \pm 10, Vacuum 20 ltr / min.
Pedal : 150 ml per stroke volume
- Pump Type** : Oil free Piston Pump (Electrical) & Piston pump (Pedal)
- Jar** : Wide mouthed autoclavable 2 ltr collection jar (Polycarbonate), safety jar with mechanical / electronic overflow motor shut off system and autoclavable bacterial filter.
- Tubing** : 8 mm ID x 2 mtr. (Silicon)
- Vacuum Gauge** : Bourden type 5 cm Dia
- Power** : Electrical: Operates on 220 V AC and Pedal (manual). In-built maintenance-free 6 V, 4.5 Amp rechargeable battery. 55 Watts.
- Dimension & Weight** : 44 x 28 x 27 cms, 9.5 Kg.
- Battery Backup** : 120 minutes on full charge.

Provided with AC adaptor and cigarette lighter socket cable for use in ambulance / car. Can be hung in ambulance or on wall.

Anand's Portable AC / DC / Pedal Suction Unit for Ambulances, Hospitals, Camps, Trains, Home care etc.

LARGEST RANGE OF SUCTION UNITS IN INDIA

HIGH VACUUM SUCTION UNITS

HI-VAC PLUSS (90 LTR/MIN)
HI-VAC PLUSS (60 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC 51
HI-VAC 30
EUROVAC-ELITE

PIPELINE VACUUM UNIT ON TROLLEY

PIPELINE VACUUM UNIT ON TROLLEY

MOBILE SUCTION UNIT

MB-36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP-30

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

ACCESSORIES

POLYCARBONATE JARS
YANKEUR SUCTION TUBE / HANDLE
KARMAN TYPE CANNULAE
SILICON VENTOUSE CUP
BACTERIAL FILTER

Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA
Phone : 91-11-25225225, 25229206, 42464264
Fax : 91-11-25225062
Email : sales@anandind.com
Web : www.anandind.com / www.suction.tel

Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Phone: 91-1276-267061, 268407

To view the product installation & usage video, go to www.anandind.com