



# ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE  
0434

## MANUAL OPERATED SUCTION UNITS ◀

P-9



### Anand's Non-Tiring Manual-Powered Suction Unit

**Housing:** Engineering plastic moulded parts - non corrosive, easy to dismantle, clean, disinfect and reassemble.

**Capacity:** Outright vacuum - 600 mm Hg  $\pm$ 10 with average per stroke volume of 200 ml.

**Pump Type:** Efficient piston pump for creating vacuum instantly.

**Jar:** Autoclavable 1 X 1 ltr. jar (Polycarbonate)

**Tubing:** 8 mm ID X 2 mtr. (Silicon)

**Vacuum Gauge:** Bourden type, 5 cms 0 - 760 mm Hg

**Dimension & weight:** 32 x 17 x 30 cms. 2.5 kg.

Portable pedal unit with top view manometer visible to the operator.

Bacterial filter (0.3 micron with single use / autoclavable housing) between jar and pump.

**Ease of operation** - one person can operate unaided. Also useful for field hospitals, midwives, medical centres, mines, aircraft, trains, mountaineers, on-foot paramedics etc.

## Product Range

### HIGH VACUUM SUCTION UNITS

HI-VAC PLUS (90 LTR/MIN)  
HI-VAC PLUS (30/65 LTR/MIN)  
MAGNAVAC  
HI-VAC  
HI-VAC (JR.)  
EUROVAC-ELITE  
EUROVAC PLUS

### MOBILE SUCTION UNIT

MB 36

### ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

### PORTABLE SUCTION UNITS

EUROVAC  
EUROLITE  
EURO-DEL

### AC / DC SUCTION UNITS

EUROVAC-A  
EUROVAC-B  
INSTAD

### AC / DC / PEDAL SUCTION UNIT

MULTIVAC

### MANUAL OPERATED SUCTION UNITS

P-9  
HANDHELD SUCTION UNIT

### ACCESSORIES

POLYCARBONATE JARS  
SILICON VENTEUSE CUP  
YANKEUR SUCTION TUBE  
KARMAN TYPE CANNULAE  
BACTERIAL FILTER

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.  
New Delhi-110026, INDIA  
Ph. +91-11-25225225, 25229206, 42464264  
Fax. 91-11-25225062

### Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Ph. 91-1276-267061, 268407