



ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

PIPELINE VACUUM UNIT ON TROLLEY ◀



PIPELINE VACUUM UNIT ON TROLLEY.

- Unit portable on 5 castors (2 with brakes)
- 2 way changeover level for 2 x 2.5 ltr PC jars.
- Additional overflow safety jars with reusable / autoclavable filter for each collection jar to prevent fluid from entering the pipeline.
- Fine vacuum regulator with vacuum gauge.
- Can be connected to any vacuum source.
- Autoclavable catheter holder provided.
- Carrying handle for easy movement.

Product Range

HIGH VACUUM SUCTION UNITS

HI-VAC PLUSS (90 LTR/MIN)
HI-VAC PLUSS (30/65 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC (JR.)
EUROVAC-ELITE
EUROVAC PLUS

MOBILE SUCTION UNIT

MB 36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

PIPELINE VACUUM UNIT ON TROLLEY

PIPELINE VACUUM UNIT ON TROLLEY

ACCESSORIES

POLYCARBONATE JARS
SILICON VENTOUSE CUP
YANKEUR SUCTION TUBE
KARMAN TYPE CANNULAE
BACTERIAL FILTER

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.
New Delhi-110026, INDIA
Ph. +91-11-25225225, 25229206, 42464264
Fax. 91-11-25225062

Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Ph. 91-1276-267061, 268407