



ANAND MEDICAIDS PRIVATE LIMITED

An ISO 13485:2003 & 9001:2000 Certified Company

CE
0434

ACCESSORIES ◀

YANKEUR SUCTION TUBE



Transparent moulded Yankeur suction handle attached with two & half meter long P.V.C. tubing suitable for easy attachment to any type of suction apparatus. A port is provided at the distal end of the handle for finger tip / Thumb control facility during suction management. Yankeur handle is open with four lateral eyes to withstand negative pressure and ensure continuous and easy suction available with / without filter.

E.T.O. Sterilized and Pyrogen free.

Product Range

HIGH VACUUM SUCTION UNITS

HI-VAC PLUSS (90 LTR/MIN)
HI-VAC PLUSS (30/65 LTR/MIN)
MAGNAVAC
HI-VAC
HI-VAC (JR.)
EUROVAC-ELITE
EUROVAC PLUS

MOBILE SUCTION UNIT

MB 36

ELECTRIC CUM MANUAL OPERATED SUCTION UNIT

EP 90

PORTABLE SUCTION UNITS

EUROVAC
EUROLITE
EURO-DEL

AC / DC SUCTION UNITS

EUROVAC-A
EUROVAC-B
INSTAD

AC / DC / PEDAL SUCTION UNIT

MULTIVAC

MANUAL OPERATED SUCTION UNITS

P-9
HANDHELD SUCTION UNIT

ACCESSORIES

POLYCARBONATE JARS
SILICON VENTEUSE CUP
YANKEUR SUCTION TUBE
KARMAN TYPE CANNULAE
BACTERIAL FILTER

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office.

No.16, Road No.33, Punjabi Bagh Ext.
New Delhi-110026, INDIA
Ph. +91-11-25225225, 25229206, 42464264
Fax. 91-11-25225062

Works.

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Ph. 91-1276-267061, 268407