



## LINER

1 Ltr., 2 Ltr. & 3 Ltr.

- Soft bag material : LDPE (low density polyethylene)
- Individually packed in blisters



## CANISTER

1 Ltr., 2 Ltr. & 3 Ltr.

- Soft bag material : LDPE (low density polyethylene)
- Individually packed in blisters

|                                 | FLOVAC®                                                |
|---------------------------------|--------------------------------------------------------|
| Lid material                    | HDPE (high density polyethylene)                       |
| VACUUM port                     | Conic connector, female                                |
| PATIENT port                    | Ø 14.0 ÷ 15.5 mm (Ø 8.0 ÷ 9.2 mm with elbow connector) |
| TANDEM port                     | Ø 8.0 ÷ 9.2 mm                                         |
| Maximum suction value           | -750 mbar (-570 mm Hg)                                 |
| Maximum flow value              | 42 LPM under recommended conditions                    |
| Maximum graduation interval     | 50 ml                                                  |
| Patient hose sizes              | inner Ø ≥ 6 mm - L max =2.5 m                          |
| Vacuum hose sizes               | inner Ø ≥ 6 mm - L max =1.8 m                          |
|                                 | <b>FILTER</b> main characteristics                     |
| Filtration efficiency (typical) | >99.999995% with particle size of 0.1µ                 |
| Membrane                        | 100% expanded PTFE GORE™ Medical Membrane              |
| Prefilter                       | micro fiberglass                                       |
| Support                         | non-woven PE/PES                                       |



Support ring for FLOVAC®

## SUCTION UNIT JARS



0.5 ltr jar  
with hook

1 ltr  
Oval shaped  
jar with hook

1 ltr jar  
with Hook

1 ltr jar

2 ltr jar  
Also Available with Hook

2.5 ltr jar  
with handle

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
Phone : 91-11-25225225, 25229206  
Fax : 91-11-25225062  
Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Phone: 91-1276-267061

[www.anandind.com](http://www.anandind.com)