

# DISEASE TRANSMISSION

Avoid it with

**Anand's** INFECTION CONTROL SYSTEM **NEED OF THE HOUR**

- Prevents contamination : Hygienic / convenient
- Maximum protection for patients, hospital staff and the suction source



## CENTRAL PIPELINE SUCTION TROLLEY

All the necessary devices included in one suction system for a broad range of surgical and medical procedures

## CENTRAL PIPELINE SUCTION TROLLEY

- Mobile Suction Unit compatible with Central Pipeline Vacuum outlets.
- Provided with compatible adapter, 4 collection jars with disposable liners (2 Ltr. each), regulator etc.
- With Bracket System for collection jar support
- Dimension: 580 x 580 x 820 mm
- Weight : Approx. 5 Kg



**Top View**

Four Interconnected Jars (2 Ltr. each with Disposable Liners) are locked with clamps



**Vacuum Gauge**  
with release valve



**Connector Pipe**

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
Phone : 91-11-25225225, 25229206  
Fax : 91-11-25225062  
Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
Phone: 91-1276-267061