


EP-30 Back View (Manual Suction)



EP-30

Main Applications

- Minor OT's
- Intensive Care Unit
- Emergency Department
- Wards



TECHNICAL FEATURES

Model	EP-30 ELECTRIC CUM PEDAL SUCTION UNIT
Housing	MS Powder Coated, SS Top on two wheels
Max. Vacuum / LPM	Electric 600 mm Hg, 45 Ltrs / min / Pedal 600mm Hg, 200 ml /push stroke
Pump Type	Electric - Diaphragm / Pedal - Piston
Jars	2 x 2 Ltr. PC Jar
Filter	Bacterial Filter Autoclavable / Reusable
Tubing	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	6.25 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz 90 Watt
Noise Level	55 dB A \pm 3
Dimensions	42 x 39 x 64 cms
Net Weight	15.5 Kg
Overflow Safety	Dual Mechanical
Optional	110 V
Certifications	CE, IEC 60601-1-2012