



## EUROLITE

### Main Applications

- Nursing / Senior Citizen
- Home Care
- Neonatal / Paediatrics

## TECHNICAL FEATURES

Model	EUROLITE	EUROLITE (FOR VACUUM DELIVERY)
Housing	ABS Plastic Moulded	ABS Plastic Moulded
Max. Vacuum / LPM	-600 mm Hg, 21 Ltrs / min	-600 mm Hg, 21 Ltrs / min
Pump Type	Piston pump	Piston pump
Jars	1 Ltr. PC Jar	1 Ltr. PC Jar
Filter	Bacterial Filter Autoclavable / Reusable	Bacterial Filter Autoclavable / Reusable
Tubing	8 mm ID x 2 mtr. (PVC / Silicon)	8 mm ID x 2 mtr. (PVC / Silicon)
Vacuum Gauge	5 cm dia 0-760 mm Hg	5 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz 90 Watt	220 V AC, 50 / 60 Hz 90 Watt
Noise Level	55 dB A ± 3	55 dB A ± 3
Dimensions	30 x 14 x 33 cms	30 x 14 x 33 cms
Net Weight	4.5 Kg	4.5 Kg
Overflow Safety	Mechanical / Electronic Acoustic Alarm	Mechanical / Electronic Acoustic Alarm
Accessories / Optional	110 V on request	Ventouse Cup / 110 V on request
Certifications	CE, IEC 60601-1	CE, IEC 60601-1



## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
 Phone : 91-11-25225225, 25229206  
 Fax : 91-11-25225062  
 Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
 Phone: 91-1276-267061