



## FLOVAC 50

### Main Applications

- Emergency Department
- Intensive Care Unit
- OT / Surgery
- Plastic Surgery
- Obstetrics
- Gastroenterology

## TECHNICAL FEATURES

Model	FLOVAC 50
Housing	Plastic moulded on two wheels
Max. Vacuum / LPM	-700 mm Hg, 60 Ltrs / min
Pump Type	Diaphragm pump
Jars	2 x 4 ltr. wide mouthed (Polysulfone) graduated jars with changeover lever
Filter	Bacterial Filter Autoclavable / Reusable
Tubing	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	6.25 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz, 580 Watt
Noise Level	60 dB A $\pm$ 3
Dimensions	46 x 34 x 86 cms
Net Weight	18 Kg
Overflow Safety	Mechanical / Electronic Acoustic Alarm
Accessories / Optional	Foot Switch / 110 V on request
Certifications	CE, IEC 60601-1

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
 Phone : 91-11-25225225, 25229206  
 Fax : 91-11-25225062  
 Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
 Phone: 91-1276-267061