

THE NEW **VACUUM DELIVERY**

SUCTION UNIT FROM **Anand**[®]

VERY USEFUL FOR THE PROCESS OF EXTRA GENTLE VACUUM ASSISTED DELIVERY

INSTAD-DEL

→ Vacuum Delivery Suction

The Foot Switch operation allows the doctor to increase the pressure steadily without touching the machine.

The required pressure (600 mm Hg) can be attained in seconds and the pressure will be maintained.



The Vacuum pressure does not drop due to a silicon non-return valve and the doctor can concentrate only on delivery





Height Trolley for easy operation



Digital Display of Vacuum

Can be adjusted from 0-600 mm Hg



Foot Switch

The Foot Switch helps the operator to increase the vacuum to required limit and to release & adjust the vacuum level from the Foot Switch

A Vacuum Buildup

M Vacuum Release



INSTAD-DEL



Low Power Consumption



Oil-Free Piston Pump



Low Noise Level 55 dba



Overflow protection in Jar



Available with **3 cups (Silicon)**
50 mm / 60 mm / 70 mm
suitable for Occipito Anterior Positions

3 Bird type cups (SS)
40 mm / 50 mm / 60 mm for Occipito Anterior

1 Bird Cup
50 mm for Occipito Posterior



ANAND MEDICAIDS PRIVATE LIMITED

No.16, Road No.33, Punjabi Bagh Ext., New Delhi-110026, INDIA
Phone : 91-11-25225225, 25229206, 25225062
Mobile : 91-9310016965
Email : sales@anandind.com

Works:
1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
Mobile : 91-9873111104