

SUCTION UNITS AT A GLANCE

	HIGH VACUUM							ELECTRICAL & MANUAL
Model	HI VAC PLUS VARIABLE	HI VAC PLUS 90 Ltr.	HI VAC PLUS 60 Ltr.	HI VAC MS / SS CABINET	HI VAC FOR LIPO SURGERY	HI VAC JR.	EUROVAC ON TROLLEY	EP-30
Housing	Plastic moulded on two wheels with foot switch & changeover lever	Plastic moulded on two wheels with foot switch & changeover lever	Plastic moulded on two wheels with foot switch & changeover lever	MS / Stainless Steel on four castors	Stainless Steel on four castors	MS / Stainless Steel on four castors	ABS Plastic Moulded on four castors with foot switch	MS Powder Coated, SS Top on two wheels
Max Vacuum	700 mm Hg	700 mm Hg	600 mm Hg	700 mm Hg	700 mm Hg	650 mm Hg	575 mm Hg	Electric -600 mm Hg Pedal -600 mm Hg
Max Flow	Variable Flow Rate 40 / 50 / 60 Ltrs / min with easy-touch button	90 Ltrs / min	60 Ltrs / min	45 Ltrs / min	70 Ltrs / min	30 Ltrs / min	45 Ltrs / min	Electric 45 Ltrs / min Pedal 200 ml / push stroke
Battery Backup /operative Time	-	-	-	-	-	-	-	-
Standard Operation (Voltage)	220 V AC 50 / 60 Hz 360 Watt	220 V AC 50 / 60 Hz 360 Watt	220 V AC 50 / 60 Hz 200 Watt	220 V AC 50 / 60 Hz 380 Watt	220 V AC 50 / 60 Hz 380 Watt	220 V AC 50 / 60 Hz 180 Watt	220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 90 Watt
Jars	2	2	2	2	2	2	2	2
Net Weight	20.5 Kg	20.5 Kg	17 Kg	34 Kg	34.5 Kg	23.5 Kg	8 Kg + 9 Kg Trolley	15.5 Kg



PORTABLE		BATTERY OPERATED		BATTERY & MANUAL		MANUAL OPERATED	
EUROVAC AC	INSTAD AC	EUROVAC AC / DC	INSTAD AC / DC	EUROVAC A	MULTIVAC	P-9	HANDHELD
ABS Plastic Moulded	ABS Plastic Moulded	ABS Plastic Moulded	ABS Plastic Moulded	MS Powder Coated	MS Powder Coated	Plastic Moulded	Plastic Moulded
600 mm Hg	700 mm Hg	600 mm Hg	600 mm Hg	600 mm Hg	600 mm Hg	600 mm Hg	450 mm Hg
30 Ltrs / min	30 Ltrs / min	55 Ltrs / min	21 Ltrs / min	21 Ltrs / min	21 Ltrs / min Pedal 200 ml / push stroke	200 ml / push stroke	-
-	-	180 minutes (with fully charged /new battery)	90 minutes (with fully charged /new battery)	90 minutes (with fully charged /new battery)	90 minutes (with fully charged /new battery)	-	-
220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 70 Watt 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC, 50 / 60 Hz, 70 Watt, 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC 50 / 60 Hz 70 Watt 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC 50 / 60 Hz 70 Watt Electric - Pedal 12 V DC Internal Rechargeable Battery 12 V 4.5 Amp	-	-
1	1	1	1	1	1	1	1
8 Kg	3 Kg	10.5 Kg	5.5 Kg	5.5 Kg	8.5 Kg	3 Kg	192 gm