



## EUROVAC-A

### Main Applications

- Hospital Crash Carts
- Patient Transfer
- Long Term Home Care
- Tracheostomy
- EMS / Ambulance
- Bed Side

## MULTIVAC

### Main Applications

- Hospital Crash Carts
- Patient Transfer
- Long Term Home Care
- Tracheostomy
- EMS / Ambulance
- Bed Side



## TECHNICAL FEATURES

Model	EUROVAC - A	MULTIVAC
<b>Housing</b>	MS Powder Coated	MS Powder Coated
<b>Max. Vacuum / LPM</b>	-600 mm Hg, 21 Ltrs / min	Electrical : -600 mm Hg, 21 Ltrs / min (30 Ltrs / min on request) Pedal : 200 ml per stroke
<b>Pump Type</b>	Piston pump	Piston pump
<b>Jars</b>	1 x 0.5 Ltr. PC Jar	1 x 1 Ltr. PC Jar
<b>Filter</b>	Bacterial Filter Autoclavable / Reusable	Bacterial Filter Autoclavable / Reusable
<b>Tubing</b>	8 mm ID x 2 mtr. (Silicon)	8 mm ID x 2 mtr. (Silicon)
<b>Vacuum Gauge</b>	5 cm dia 0-760 mm Hg	Bourdon type 5 cm dia 0-760 mm Hg
<b>Power Supply</b>	220 V AC, 50 / 60 Hz, 55 Watt, 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC, 50 / 60 Hz, 55 Watt, Electric - Pedal 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)
<b>Noise Level</b>	55 dB A $\pm$ 3	55 dB A $\pm$ 3
<b>Battery</b>	Lead Acid (Ni-cd Optional)	Lead Acid (Ni-cd Optional)
<b>Battery Backup</b>		120 minutes (240 minutes on request)
<b>Dimensions</b>	37 x 14 x 30 cms	33 x 30 x 41 cms
<b>Net Weight</b>	5.5 Kg	8.5 Kg
<b>Overflow Safety</b>	Mechanical / Electronic Acoustic Alarm	Mechanical / Electronic Acoustic Alarm
<b>Accessories / Optional</b>	Car Charger & Power Adaptor (Included)	Car Charger & Power Adaptor (Included)
<b>Certifications</b>	CE, IEC 60601-1	CE, IEC 60601-1

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
 Phone : 91-11-25225225, 25229206  
 Fax : 91-11-25225062  
 Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
 Phone: 91-1276-267061

[www.anandind.com](http://www.anandind.com)