



India's Most Trusted Suction Brand




VARIABLE FLOW RATE



40 / 50 / 60 LPM with easy-touch button

CHANGEOVER LEVER



ELECTRONIC SENSOR



JAR



HI VAC PLUS VARIABLE

- Main Applications**
- Emergency Department
 - Intensive Care Unit
 - OT / Surgery
 - Plastic Surgery
 - Obstetrics
 - Gastroenterology

TECHNICAL FEATURES

Model	HI VAC PLUS - VARIABLE
Housing	Plastic moulded on two wheels & changeover lever
Max. Vacuum / LPM	700 mm Hg 40 / 50 / 60 LPM with easy-touch button
Pump Type	Oil free Piston pump
Jars	2 x 5 Ltrs. PC Jar
Filter	Bacterial Filter Autoclavable / Reusable
Tubing	10 mm ID x 2 mtr. (PVC)
Vacuum Gauge	6.25 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz, 360 W
Noise Level	55 dBA ± 3
Dimensions	46 x 34 x 86 cms
Net Weight	20.5 Kg
Overflow Safety	Mechanical & Electronic Sensor
Optional	4 Ltr. Polysulfone (PSU) Jars / 110 V / Foot Switch
Certifications	CE, IEC 60601-1-2012