



India's Most Trusted Suction Brand



HANDHELD

Main Applications

- Patient Transfer
- EMS / Ambulance
- Home Care



P-9

Main Applications

- Patient Transfer
- Emergency Department
- Home Care

TECHNICAL FEATURES

Model	P-9 MANUAL SUCTION UNIT	HAND HELD MANUAL SUCTION UNIT
Housing	Plastic Moulded	Plastic Moulded
Max. Vacuum / LPM	600 mm Hg, 200 ml / push stroke	450 mm Hg
Pump Type	Piston Pump	-
Jars	1 x 1 Ltr. PC Jar	1 x 200 ml Jar
Filter	Bacterial Filter Autoclavable / Reusable	-
Tubing	8 mm ID x 2 mtr. (Silicon)	-
Vacuum Gauge	5 cm dia 0-760 mm Hg	-
Power Supply	N / A	N / A
Noise Level	N / A	N / A
Dimensions	33 x 18 x 30 cms	17 x 5 x 8 cms
Net Weight	3 Kg	192 gm
Overflow Safety	Mechanical	N / A
Accessories	Tool Kit, Filter Paper & Push-in Nozzle	Tube
Certifications	CE, IEC 60601-1-2012	CE, IEC 60601-1-2012